



Endura-Seal Injection Port Septa



<u>Description</u>	<u>Item No.</u>
Endura-Seals Septa – 11 mm	250111
Endura-Seals Septa – Shimadzu Style	250100

Warnings

Warnings must be read carefully and understood. Improper use of this product can cause harm or death to personnel and damage to property!



Wear eye protection and use caution when working with pressurized systems.



When using a high-pressure gas source overpressure protection must be provided.

Injection Port and the Endura-Seal nut may be hot after use.

Max Pressure: 100 PSI. Max Temp 400C.

For use with 23 gauge cone tip needles only.

Special precautions might be required when using hydrogen or hazardous gases. Consult local regulations and your company's safety procedures.

Setup and Use

Endura-Seals Installation

1. Start with a clean injection port, and inlet pressure at zero and turned off.
2. Gently press the Endura-Seal into the injection port.
3. Screw on the included nut to finger tight plus ¼ turn. Do not overtighten.



Trajan Scientific Americas Inc
2601 Technology Drive,
Louisville, KY 40299, USA
Tel: 502 491 6300
Fax: 502 491 3390
usa@trajanscimed.com



Recycling

Please contact CRS for information regarding recycling and reconditioning gas purifier products.

© Trajan Scientific Americas Inc
09/2022 995285 _RevB