

Accessories | Manual

Our ergonomically designed crimping tools have advanced the standards of manual crimpers and decappers by addressing all of the deficiencies associated with traditional blocky, metal tools.

- Reliably seal crimp caps manually with less physical strain.

Features and benefits

- Constructed of tough, light-weight plastic that is incredibly durable, but weighs much less than metal crimping tools.
- The curved handles create a natural grip, greatly improving comfort over metal grips.
- Unique bottom-pull technology allows for steadier operation in comparison to top-push handle designs.
- The adjustment knob is easily viewed and clearly displays + and - symbols along with directional arrows, simplifying the process of adjusting the crimp setting - no tools needed!
- No more extra squeeze required! When the adjustment knob touches the body of the crimper, the desired crimp setting has been reached - significantly reducing arm strain and fatigue associated with traditional manual crimping.

Recommended applications

Low volume crimping/decapping applications for aluminum crimp caps.

Product specifications

Tools in 8, 11, 13, and 20 mm crimping and decapping.

For more information about this product visit www.trajanscimed.com or contact techsupport@trajanscimed.com.

